

Belmont Handpieces Catalog

Belmont Handpieces
Catalog

Belmont handpieces - your complete dental solution:

With the Japan-made:

- Enhance your experiences at work
- Confirm your confidence in quality
- Improve reliability on the results



For the beauty of dentistry.

Symbol



Handpiece, air motor
with illumination



Coupling with
LED illumination



This product is autoclavable
up to max.135°C.



This product can be washed
via Thermo Disinfectant.



Clean Head System

Air Turbines

N-5 Series
DA-391 Turbine

• Stainless Steel Body



N-5 Series
DA-312 Turbine

• Stainless Steel Body




N-2 Series
DA-520 Turbine

• Stainless Steel Body



DA-385 LED Coupling

• Titanium Coating • Stainless Steel Body
• Non-Retraction Valve



DA-525 Coupling

• Stainless Steel Body • Non-Retraction Valve



DA-526 Coupling

• Non-Retraction Valve



Air Motors

N-5 Series
DA-281 Air Motor

• Torque 1,5 N-cm and more



N-5 Series
DA-290C Contra-Angle

• 1:1 Direct Drive



N-5 Series
DA-290S Straight Handpiece

• 1:1 Direct Drive



Category	Grade	Model	Type	Optic	Power	Spray	Speed (min ⁻¹)	Head Size	Bur Type	Chuck	Bearing	Weight [g]
Air Turbine	N-5	DA-391	Standard Head	LIGHT	20W	4x	360,000-430,000	Ø12.1 x H13.3mm	ø1.59-1.60mm Std	Push Button	Ceramic	Approx. 46
		DA-312	Standard Head	NO LIGHT	20W	4x	360,000-430,000	Ø12.1 x H13.3mm	ø1.59-1.60mm Std	Push Button	Ceramic	Approx 45
	N-2	DA-520	Standard Head	NO LIGHT	17W	4x	350,000-430,000	Ø12.1 x H13.3mm	ø1.59-1.60mm Std	Push Button	Ceramic	Approx 51
Air Motor	N-5	DA-281	Internal Water Spray	NO LIGHT	-	-	22,000±10%	-	-	-	-	Approx 64
		DA-290C	1:1 Contra-Angle	NO LIGHT	-	1x	Max 40,000	Ø9.9 x H13.8mm	Ø2.35 CA	Push Button	-	Approx 80
		DA-290S	1:1 Straight Handpiece	NO LIGHT	-	1x	Max 40,000	-	Ø2.35 HP (CA*)	-	-	Approx 81

* When using the bur stopper [Option].
The specification and colors are subject to change without notice.